WOOKVIBE.COM

**EN** Instructions manual  
**SI** Navodila za uporabo  
**CRO** Upute za uporabu

WOOKVIBE ELITE | WOOKVIBE PLUS | WOOKROLL | WOOKIE

To reduce the risk of malfunction, electric shock, injury, fire, death, and personal injury - read all of the instructions and warnings in this manual carefully before using the devices in these instructions.

## Precautions!

1. Consult your doctor before using the devices in these instructions or in case of muscle injury.
2. It is forbidden to use the devices on the surface of non-dry and non-clean skin. The devices can be used over clothing where you gently press and move over the muscle.
3. It is forbidden to use the devices on parts of the body, above the neck and on parts with less muscle coverage, or directly on the bones. Use them only on the soft tissues of the body without pain or discomfort.
4. Do not use the devices if bruising, pain, or discomfort occurs.
5. The devices must not come into contact with water.
6. Do not place appliances on fire or near other heat sources. Store devices in a dry place.
7. It is forbidden to intentionally damage devices.
8. Children between the ages of 6 and 18 must use the devices under the supervision of parents or guardians.
9. Do not disassemble the devices in any way.
10. Do not hold your fingers, hair, or other parts of your body near massage attachments or holes, as you may injure yourself.
11. After charging, unplug the charger.
12. Use the devices according to the instructions.
13. We recommend the included chargers / cables to charge the battery.
14. Inspect your device carefully before each use.

## Additional information

Read the instructions manual carefully before using the devices and use the devices correctly according to the instructions. Our company will not be liable for any damage or loss to the customer resulting from improper use of the device as explained. Do not use the devices in any of the following situations without the prior consent of your doctor:

- complications of pregnancy and urination (such as neuropathy or retinal damage),
- wearing a pacemaker,
- recent operation,
- epilepsy or migraine,
- intervertebral disc herniation,
- spondylolisthesis,
- vertebroarterial deviation or solid joint arch,
- recent joint replacement operations,
- metal needles or steel plates associated with your health.

These contraindications do not mean that you can not use this device. We suggest that you consult your doctor specialist or physiotherapist before use.

## **Safety instructions:**

1. Use the massage gun only as directed. Read and follow the instructions that came only with your purchase. Do not open or repair the massage gun and use only the attachments provided.
2. Devices are not suitable for children, for people who have problems with reduced sensory perception and who suffer from mental health problems. The devices are not toys, so they are not intended for children to play with. If you use it near children, you need to be very careful.
3. When the devices are in operation, do not attempt to remove or tamper with the battery.
4. Use only the supplied charger. A different charger may cause personal injury or damage. Never plug the devices into an electrical outlet or any outlet in the car. This can permanently damage the battery.
5. Charge the devices indoors in a well-ventilated and dry place. Do not charge them outside, in the bathroom or near a bathtub or pool. Do not use the devices or chargers in wet areas. The charger should never be exposed to water, moisture, rain or snow.
6. Never use the devices near explosive objects (gases, dust or flammable materials) as there is a risk of sparks causing a fire.
7. Ensure proper handling of your charger, do not leave it always on. Unplug the charger from the electrical plug and the device when not in use for a long time. When charging the devices, be careful how you leave the cord so as not to obstruct anyone. If the charger or cable is damaged, it must be replaced immediately.
8. Do not block any ventilation on your device and do not use it if any aerator is blocked. The devices should be free of dust, hair, fibers or other objects that can reduce air flow. Do not insert hair, pieces of clothing and fingers into the aerators and moving parts of the devices.
9. Never attempt to burn or expose the devices or battery to a fire. This will cause personal injury and trigger toxic fumes from burning the battery.
10. Don't overdo the charging. It is recommended to stop charging when the devices are fully charged. Overcharging will adversely affect battery life.
11. Do not damage the charger. If your charger has been hit hard, dropped, or damaged, do not use it.
12. Never use your devices under any bedding, blanket or pillow. This can cause overheating and can cause a fire.

## **Charging instructions**

1. Fully charge your device for a few hours before using it for the first time.
2. Use the supplied charger when charging.
3. The LEDs on the device flash regularly, indicating battery power and indicating charging.
5. The battery can be charged at any time and at any level of battery. Do not use the device when charging.
6. When the battery is at less than 25% battery level, the device needs to be charged.
7. The average operating time of a product with a full battery is about 5-6 hours, depending on the level of use.

## Instructions for use of WookVibe massage guns

1. Select the desired massage attachment, insert it into the massage gun and turn on WookVibe. Do not place the massage gun on your body. You may find it as an insignificant, small detail, but it is important for better control of the WookVibe massage gun.
2. Before embarking on the desired therapy, remember to always start at lower power and gradually increase the level of power. Too much power can lead to additional muscle tension.
3. When you get to a knot, stiffness, or muscle pain site, do not massage immediately directly at that site, but start a few inches away and then slowly move toward the aforementioned point.
4. Do not press immediately, but let the WookVibe massage gun slide over the muscle. Repeat this a few times and then gradually add a little pressure or increase the power of the massage gun.
5. When you feel that the muscle has already adjusted to the power level, leave the massage gun at the point of knot, pain or stiffness for 1-2 minutes and relax.

Then repeat the above steps for the other points on the body as well.

### Additional information:

- Massage each muscle for a maximum of 3 minutes.
- For a more effective massage, replace the massage attachment suitable for the muscle group
- Breathe deeply to relax your body and muscles.
- Listen to your body - when in pain, reduce the pressure or adjust the power level

## Can I use the WookVibe massage gun on all parts of the body?

A massage gun can be used on the shoulders, buttocks, back, neck, thighs and hips. It works by massaging your body tissues and thus relieving muscle pain and smoothing muscle knots that could cause muscle pain and stiffness. Do not use the massage gun on the face or on the head area.

There are various injuries and conditions in which the use of a massage gun is not recommended. For certain muscle strains and injuries, a massage gun will not help and can do more harm than good. If you are not sure if you need to use it for a specific injury, consult your doctor or physiotherapist.

Injuries in which you must avoid using a massage gun:

- Pulled muscles - Do not use the massage gun on already injured, pulled muscles. Massage and stretching can only further damage such muscles.
- Stretched ligaments - This happens when you stretch a muscle outside of its normal range of motion. When this happens, you usually hear a bang-like sound.
- Inflammations such as tendonitis, bursitis, fasciitis, periostitis, etc.
- Bone fractures - It may seem obvious, but don't use a massage gun anywhere near or on broken bones. Even if you no longer have pain, the area is still weak, so always consult your doctor or physiotherapist before using a massage gun.

## Adjusting the massage attachment

Hold the massage gun in your hand, align the nozzle on the mouth of the massage gun, slowly insert the nozzle (direction in Figure 1) and position it as shown in Figure 2.

Figure 1: direction of insertion of massage attachments

Figure 2: Properly inserted massage attachment

## How long should I massage a certain area?

It depends on the needs of a particular muscle group. Once you decide how long to use a massage gun in a particular area, you need to use your own judgment. While using a massage gun, you will encounter stubborn muscle knots that need extra time to ease.

The recommendation is a minimum of 15 seconds and a maximum of 3 minutes on any muscle group. You may want to hold the massage gun at a certain point for a longer period of time if you feel the need.

**30 second muscle activation** - use a massage gun before you start exercising or in the morning. Awaken the muscles to be ready to exercise.

**15 second muscle reactivation** - during exercise, where you use a massage gun for 15 seconds on a selected muscle group. This will reactivate the muscle so that you get tired more slowly and last longer under exertion.

**2 to 3 minutes of pain relief and regeneration** - 2 to 3 minutes after the workout, massage with a massage gun over any inflamed, tired, painful or stiff muscle group.

## Maintenance, cleaning, storage and transport

1. Clean the device with a plain or paper towel.
2. In case of storage or travel, place it in the enclosed carry case.
3. It must be stored in a dry and well-ventilated place free of corrosive gas.
4. If not used for a long time, the device shall be cleaned, stored and charged every six months.
5. The product must be secured during transport to prevent vibrations, rolling and falls which could potentially damage the product.

## Massage gun WookVibe Elite


1. Powerful brushless motor with quiet sliding technology.
2. Metal heat sink.
3. Interchangeable massage attachment.
4. Battery. The battery is not removable. Do not remove it by force. Charge the battery before first use.
5. LCD screen with adjustment button. Nine levels of intensity and two massage modes.
6. Silica gel handle that is smooth and comfortable to use.
7. Charging hole. Charge the product with the supplied charger.

Dimensions: 21 x 16.5 x 4.5 cm;  
 Weight: 0.80 kg  
 Color: black, gray  
 Certificates: CE, FCC, UN38, EMC, PSE, MSDS;  
 Measured el. power: 60W;  
 El. voltage: 14.4-16.8V;  
 Battery: 2600 mAh Lithium-ion;  
 Speed (rounds per minute): 1200-3500;  
 Usage time: up to 360 min (7 days of use per 1 person);  
 Charging time: up to 120 min;  
 Massage attachments: 6  
 Warranty: 1 year

1. "FITNESS" mode: High speed and high torque. With 9 ascending levels.
2. "MEDI" mode: pulsating power and speed, suitable for therapy purposes
3. Power button: hold for two seconds to turn on and off, short press to shift gear. Press twice quickly to switch between FITNESS and MEDI modes.
4. Level indicator: 1-9 levels.
5. Battery status indicator. When the battery status light is blinking, it indicates that the battery status will soon be reduced.

## Massage gun WookVibe Plus

1. Battery status indicator. When the light is on at 100%, the battery is charged. When the light is on at 80%, 60%, 40%, 20%, the battery power gradually decreases. When the "Power" light is on, the battery will soon discharge and need to be charged.

2. Level indicator: on the left side of the screen, the number 20 Hz represents the 1st level, the frequency of which is 1200 rpm. If you increase the speed, the frequency will increase by 575 rpm. The highest is the 6th level - 3500 rpm.

3. Power button.

4. +, - buttons at the bottom: to increase or decrease the power level. If you press the + button, the level will increase by one, to -, the level will decrease by one level.

1. Powerful brushless motor with quiet sliding technology.
2. Metal heat sink.
3. Interchangeable massage attachment.
4. Battery. The battery is not removable. Do not remove it by force. Charge the battery before first use.
5. Power level adjustment keys.
6. Silica gel handle that is smooth and comfortable to use.
7. Charging hole. Charge the product with the supplied charger.


Dimensions: 21 x 16.5 x 4.5 cm;  
 Weight: 0.80 kg  
 Color: white, red  
 Certificates: CE, FCC, UN38, EMC, PSE, MSDS;  
 Measured el. power: 60W;  
 El. voltage: 14.4-16.8V;  
 Battery: 2600 mAh Lithium-ion;  
 Speed (rounds per minute): 1200-3500;  
 Usage time: up to 360 min (7 days of use per 1 person);  
 Charging time: up to 120 min;  
 Massage attachments: 6  
 Warranty: 1 year

1 Battery status indicator

2 Level indicator

3 Power button

4 Buttons +, -

## Massage attachments and which muscles they are intended for

Each attachment has its own function, but this does not mean that some attachments, which are in principle intended for smaller muscle groups, should not be used on larger muscle groups and vice versa. Use attachments that you think are appropriate to achieve the goals of the necessary therapy.

For relaxation of larger muscle groups. Recommended use for massaging more sensitive muscles.

For larger muscles and to help speed up blood circulation.

To reach the smaller muscles around the joints and the deeper muscle tissues.

For the spine, shoulders and around the Achilles tendon. Use on smaller stages.

For more precise massaging of smaller muscles.


Suitable for the whole body. The most versatile massage attachment.

## Use of the WookRoll vibrating massage roller

1. Move the switch from OFF to ON and place the WookRoll vibrating roller above or below the muscle you intend to handle.
2. After settling into a comfortable position, press the center button. Blue lights will light up to indicate device operation and current power level. The green lights indicate the battery status.
3. Continue to press the center button (do not hold it) until you reach the desired vibration speed.
4. When you want to turn off the vibration, simply turn the top button to OFF and wait until the blue lights turn off or until the vibration stops.

## Charging instructions

When charging the device, use the cable supplied with the roller. Then plug the cable into any USB device or USB charger (up to 5.0V, up to 2.0 A). The roller charges when the lights are flashing and the device is full when the lights stop flashing.

## How to use the WookRoll massage roller?

Always massage on muscle tissue and do not roll over tendons, joints or bone structures. Do not roll on areas that are too painful or do not roll "smoothly". We recommend that you start by putting the roller on a sensitive location or knot and gradually increase the pressure on the roller or the vibration speed. Stay on the sensitive area until the muscle relaxes, but no longer than one minute. Always use the roller before stretching.

## Use of the Wookie massage vibrating ball

1. Press the power button and place the Wookie massage vibrating ball above or below the muscle you intend to treat.
2. After settling into a comfortable position, press the same button to increase and decrease the vibration level. Red lights will light up to indicate device operation and current vibration level.
3. Continue to press the same button (do not hold it) until you reach the desired vibration speed.
4. When you want to turn off the vibration, press the same button until the Wookie massage ball turns off.

Specifications:  
Weight: 372 g  
Charging time: up to 4 hours  
Engine speed (RPM):  
2000 - 1st level,  
4000 - 2nd level,  
6000 - 3rd level  
Pulsating - 4th level  
Material: ABS + silicone  
Battery: Lithium 3.7V 2500mAh  
Operating time: 60 minutes (6 cycles), 10 minutes per cycle

## Charging information

The product can be charged with various adapters that have from 3.7 to 5.0 volts. You can use a laptop or desktop computer, a mobile powerbank, or a mobile phone charger.

Charging time is 2-4 hours. The actual charging time depends on the remaining battery and at what stages it is used.

During charging, the power and standby indicator lights will flash blue and red.

The device must not be switched on during charging. When the device is full, the battery indicator light will be red.

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